


Kod: **P**

Prelungitor



| Cod   | A    | B[mm] | C[mm] | L [mm] | Ambalare individual |
|-------|------|-------|-------|--------|---------------------|
| P10C  | G1/2 | 10    | 8,5   | 10     | 10                  |
| P110C | G3/4 | 12    | 8     | 10     | 1                   |
| P15C  | G1/2 | 10    | 12    | 15     | 10                  |
| P115C | G3/4 | 12    | 13    | 15     | 1                   |
| P20C  | G1/2 | 10    | 12    | 20     | 10                  |
| P120C | G3/4 | 12    | 15    | 20     | 1                   |
| P25C  | G1/2 | 10    | 12    | 25     | 10                  |
| P125C | G3/4 | 12    | 15    | 25     | 1                   |
| P30C  | G1/2 | 10    | 12    | 30     | 10                  |
| P130C | G3/4 | 12    | 15    | 30     | 1                   |
| P40C  | G1/2 | 10    | 12    | 40     | 10                  |
| P140C | G3/4 | 12    | 15    | 40     | 1                   |
| P50C  | G1/2 | 10    | 12    | 50     | 10                  |
| P150C | G3/4 | 12    | 15    | 50     | 1                   |